


InBody

See what you're made of


What is Body Composition?



HEIGHT
06 00
FEET INCHES


WEIGHT
250.0
POUNDS

PBF **PBF**
31% 10%




Muscle vs. Fat Density

1 lb. of muscles takes up less space than 1 lb. of fat





1lb. of fat



1lb. of muscle

Revolutionizing BIA Technology with InBody

InBody's medical-grade body composition analyzers rely on four pillars of technology to give you extremely accurate and precise BIA results that are highly correlated to gold-standard methods.




Not a product,
But a **Solution**





InBody Dial

Tracking your progress with a normal weight scale is limited. What really matters is body composition. Increasing muscle and losing fat makes your body slimmer, metabolically active and healthy. The InBody Dial changes the paradigm of a normal weight scale, starting from percent body fat, muscle mass, BMI to visceral fat level.



PBF

Percent Body Fat



SMM

Skeletal Muscle Mass



BMR

Basal Metabolic Rate



VFL

Visceral Fat Level



BMI

Body Mass Index

InBody Dial

After every test, your results are recorded and organized in the InBody app. See your test history, spot trends and share your results with your personal trainer or doctor. Make adjustments based on accurate data, so that you are consistently trending toward the right direction.



ACCURATE AND PRECISE RESULTS

Get results you can trust with the smart weight scale created by the founder of modern BIA technology.





DATA THAT MATTERS

You'll get both your Percent Body Fat and BMI values so you can see how PBF is more accurate and reliable than BMI.




CONVENIENT TRACKING OVER TIME

Automatically view your test results on your phone in real-time. Check out your progress anywhere, anytime.



InBody Dial



Handle/Electrode

User holds the device with hands and fingers.

Display Screen

Height, weight, percent body fat, muscle mass, and Visceral fat level is displayed on the screen.

Dial Button

Use to power On/Off, input height and check results..

Foot Plate

User stands on the foot electrode to measure lower body composition.




10:45




InBody

Results History Ranking

< 03.16.2019 (Sat) 10:55 >



Obesity Analysis



Research Parameters

Basal Metabolic Rate 1316 kcal







2:42



InBody

Results History Ranking

01.14.19 ~ Recent Test Date / All Unit Data




2:45



InBody

Results History Ranking


InBody Score



Friend Ranking	Points
1 상화님(인바디) 3day(s) ago	80
2 Me 118day(s) ago	78

Friends List Edit





10 Second Tests

It's as simple as standing on the device like a regular weight scale and holding the handlebar. Advanced technology with ease of use.

Accurate and Precise Results

Get results you can trust with the smart weight scale created by the company that pioneered modern BIA technology.

Data that Matters

You'll get both your Percent Body Fat and BMI values so you can see how PBF is more accurate and reliable than BMI.

Convenient Tracking Over Time

Automatically view your test results on your phone in real-time. Check out your progress anywhere, anytime.

InBody Dial Product Line



InBody Dial H20B
- Soft White



InBody Dial H20N
- Oatmeal Beige



InBody Dial Black
- Glimming Pearly Black

InBody Dial

- Home use – home fitness, home training
- Small fitness centers
- Medical use
- Day to day care
- Wellness centers
- Rehabilitation centers
- Spa and aesthetics